

The extended wedge method:

Atomic force microscope friction calibration for improved tolerance to instrument misalignments, tip-offset, and blunt probes

H.S. Khare, D.L. Burris, Review of Scientific Instruments 84 (2013) 055108 [[web](#)]

Instructions for using the LFM template

Dare to be first.



**MATERIALS TRIBOLOGY
LABORATORY**

This document provides instructions for using the solver template for the extended wedge method calculations, outlined in the *Review of Scientific Instruments* paper. Examples shown in the following pages correspond to MS Excel 2010.



Before you begin

Check version of MS Excel


Ensure that you are using the correct file for your version of MS Excel. The solver template is written as an MS Excel macro and has been tested to work with MS Excel 2010 and MS Excel 2007, on platforms running both 32bit and 64bit Windows Vista/7.

Enable macros

Ensure that macros are enabled in MS Excel. When you open the LFM template file, MS Excel may prompt you to enable macros – select ‘yes’ if prompted.

Installing solver add-on

The Solver add-on for MS Excel is typically installed automatically with the installation of MS Office. If however the Solver add-on has not been installed on your system, you will be prompted to install this when you open this file (typically in MS Excel 2007). If prompted, select ‘yes’ and follow on-screen instructions to install this add-on. Such an installation will typically not require you to insert the original MS Office CD.



Step 1: After opening the LFM Template file, you may be shown a security warning indicating that macros are disabled. Click 'Enable Content' to make sure macros are enabled.

Clipboard Font Alignment Number Styles Cells

E8

Directions: wedges 1 and 2 do not need to be equal and opposite
 prime denotes slope with respect to normal voltage
 R/h can be estimated with manufacturer reports
 Input AFM data into gray cells
 Use the solver function (data tab; add-in) to iterate output variables until residual is minimized

	inputs				outputs		
	angle (°)	Δ'	W'	R/h	μ	S	O/h
wedge 1	-6.511	-0.115	0.037	0.000	0.190	5.128	0.479
wedge 2	6.334	-0.070	0.039		0.200		
	Eq. 11	5.195					
	Eq. 12	5.195					
	Eq. 13	5.195					
	Eq. 14	5.195					
	residual	0.000					

Show dialog box

Extended Wedge Method Solver

Angle 1 (degrees)

Angle 2 (degrees)

$W'1$

$W'2$

$D'1$

$D'2$

R/h

$\mu 1$ O/h

$\mu 2$ S

Residual

Step 4: The values of friction coefficients $\mu 1$ and $\mu 2$, tip-offset ratio O/h and sensitivity S are evaluated and displayed both in the dialog box, as well as the data spreadsheet. The value of the residual (standard deviation in the values from Eqns. 11-14) is also shown.

A rectangular button with a light gray background and a thin black border. The word "Clear" is centered in the button in a black, sans-serif font.

The 'Clear' button clears the fields of the user-input dialog box

A rectangular button with a light gray background and a thin black border. The word "Close" is centered in the button in a black, sans-serif font.

The 'Close' button closes the user-input dialog box

A rectangular button with a light gray background and a thin black border. The text "Show dialog box" is centered in the button in a black, sans-serif font.

Clicking on 'Show dialog box' launches the user-input dialog box